


# EG00751

## Nine Channel Precision Temperature Module



# °C

[oem.medlab.eu](http://oem.medlab.eu)

- Nine independent channels, plus reference channel
- Small size (84.5 x 44 mm ), maximum height 10 mm
- Single power supply (5 V DC), 100 mW
- RS232 and TTL level serial asynchronous interface
- Integrated patient isolation according to IEC601-1, class CF
- Packet oriented transmission with checksum
- Transmits up to 10 measurements per second
- Resolution 0.01 °C
- Measuring range 15.00 °C up to 59.00 °C
- Compatible with all YSI 400 Probes
- Defibrillation protection
- Accuracy: +/- 0.01 °C to +/- 0.02 °C + Probe Tolerance
- Made in Germany

Applications:

Patient Monitors

Temperature measurement for incubators

Industrial precision temperature measurement



# medlab

medizinische Diagnosegeräte GmbH

Helmholtzstrasse 1

D 76297 Stutensee near Karlsruhe

Germany

Tel.: +49 (0) 7244 74110 0

FAX: +49(0) 7244 74110 28

[sales@medlab-gmbh.de](mailto:sales@medlab-gmbh.de)